

SPS

Sterile Protection Shield - Head Belt System

A lightweight and sterile faceshield against the transfer of microorganisms, body fluids and particles.



- / Lightweight
- / Economic
- / Versatile
- / Sterile



Protect your patient from harmful fluids and particles


The SPS Sterile Protection Shield Head Belt System, is a medical device worn by medical personnel during surgical procedures in the operating room. It is made to protect the patient from harmful exposure of body fluids or face particles from the operating team members. Its low weight and high level of comfort makes it a lightweight but reliable barrier between the surgeons' face and the patient against the transfer of body fluids and particles.

Comfort Pads

Easy to replace foam pads for **comfort** and **hygiene**.

Lightweight Head Belt

The patented re-usable head belt is **easy to clean** and available in sizes XS - XL and M - 4XL.



Face Shield

The sterile protection shield is **ETO sterilized**, disposable, **single-use** and easy to change. The polycarbonate lens provides a maximum field of vision without any distortions or reflexion. The face shield can be opened and closed.

Headband Adjustment Knob

Easy **one hand adjustment** for comfortable and firm fit.

/ how to buy

Code	Item
SP100	SPS Head Belt, size XS - XL
SP101	SPS Head Belt, size M - 4XL
SP250	SPS Sterile Protection Shield (24 per case)
SP774	SPS Replacement Comfort Pads (1 set)

/ where to buy



To find a THI distributor in your area, visit our website: www.thigmbh.at