

## BATTERY HOLSTER HW

Connects REF 60402 Battery Pack HW with all helmets of the ViVi® and Oxygen Aerolyte series.

### Features

- Lightweight, simple to use battery connector
- Secure mechanical fixation of the battery (no drop)
- Integrated beltclip to attach to waistband or belt
- Fits easily into your pocket.



### Intended Use:

This battery holster is a component of the ViVi® System. Worn on the body under sterile protection gowns, it is used by healthcare personnel to connect the battery with the helmet and by doing so, to ensure low voltage DC power supply to the helmet.

### Specification

Dimensions	
Length	121mm
Width	83mm
Thickness	38mm
Weight	58g
Contacts	+, D, C, T, -
Housing Material	PA-12 (Polyamide)

For use with:	
Battery	60402 Battery Pack HW
Helmets	80100 ViVi® Helmet
	80120 ViVi® Helmet w/Headlight
	70100 Oxygen Aerolyte Helmet XS-XL
	70110 Oxygen Aerolyte Helmet M-4XL
	70120 Oxygen Aerolyte LED Light Helmet XS-XL
	70130 Oxygen Aerolyte LED Light Helmet M-4XL

**THI Total Healthcare Innovation GmbH**  
 Gewerbestrasse 4, 9181 Feistritz, Austria  
 phone: +43 4228 30100  
 fax: +43 4228 30100-20  
 web: www.thigmbh.at  
 mail: office@thigmbh.at